DE • ENABOND CX

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EN • ENABOND CX

Cavity washing solution with 2% Chlorhexidine

COMPOSITION
2% Chlorhexidine digluconate, deionised water.

PROPERTIES

The product is a water solution with a base of 2% chlorhexidine digluconate formulated to clean and rehydrate cavities prepared for fillings with composite resins. The use of the solution is recommended after etching with gel. Recent research has shown that chlorhexidine removes dentine debris without interacting with the sealing capacity of the dentine adhesives, and helps eliminate post-operative sensitivity caused by residual bacteria and their endo- toxins. The use of the product for rehydrating dry dentine before bonding allows a higher bonding strength than when wetting agents are not used.

INDICATIONS

Cleansing and rehydrating dental cavities prepared for filling with composite resins.

INSTRUCTIONS FOR USE

- 1) Position the rubber dam to ensure adequate isolation of the operating field.
- 2) Prepare the cavity, then rinse with water and dry with a jet of air.
- 3) Etch the enamel and dentin with 37% phosphoric acid gel. The preparation is etched in a different way: the etching gel is applied first on the enamel for 15 sec., followed by application of the etching gel on the dentine surface for 15 sec. It is imperative for the timing of these steps to be respected precisely, and not guessed approximately. The etching gel must be washed abundantly with jet of an air-water syringe (for about 5 sec.), then dried with a jet of air, ensuring that the dentine surface is not absolutely dry, but at the same time there is not an excess of humidity. If a self-etching adhesive is used, skip point 3 and proceed to point 4.
- 4) Apply the 2% chlorhexidine solution which acts as a rehydrating agent in the event of excessive drying in the previous step. In particular, place a cotton pellet soaked in the 2% chlorhexidine solution in the prepared cavity for 30 sec. The cavity must then be dried gently with the cotton pellet or with a jet of air for 1 second, without rinsing.
- 5) Apply the primer + adhesive, or only adhesive, or self-etching adhesive, on the surface of the preparation with a microbrush, rubbing energetically for a sufficient time (see the *Manufacturer's instructions*) to ensure that the dentine surface is properly covered. Then dry well with maximum air pressure for about 10 seconds by suction or by a jet of air with a syringe, then light cure (10 sec.).
- 6) Apply the flowable composite resin on the bottom of the cavity, then layer the composite paste and light cure (for 10 sec. see the *Manufacturer's instructions*) in a series of relatively small increments, keeping the tip of the curing light as close as possible to the composite material.
- 7) Check occlusion after removing the dam, then finish with low-abrasive silicone tips and finally polish and shine with brushes and respectively diamond and aluminium oxide pastes.

RISKS

Chlorhexidine digluconate can irritate the eyes and skin.

PRECAUTIONS

Use suitable protective clothing, gloves and eyewear.

WARNINGS

If the product comes in contact with the eyes, wash with plenty of water; consult your doctor if the irritation persists. If it comes in contact with the skin, wash with soap and water. Refer to the safety data sheet for additional risks and information about tolerability.

DO NOT INJECT - DO NOT SWALLOW

Keep out of the reach of children
Do not dispose of the container in the environment after use.
Do not use the product after its expiry date

CONTRAINDICATIONS

Known individual hypersensitivity to chlorhexidine.

STORAGE

The product must be stored at ambient temperature (20-25 °C), away from direct light in a dry place.

PRESENTATION

1 multidose 10 ml bottle of 2% Ena Bond CX solution in a box.
Validity 24 months.

FIELD OF USE

Reserved for professional use in dental practice.

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Micerium S.p.A.
Italméd S.r.l.